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**CE** 2797

## Photodynamic Therapy System

### User Manual for the FLUXMEDICARE® System

MDB TEXINOV®  
56 route de Ferrossière  
38110 SAINT DIDIER DE LA TOUR  
France

Tel.: +33 (0)4 74 94 44 75  
Email: medical@texinov.fr

EN

# Contents

---

<b>1.</b>	<b>Introduction .....</b>	<b>32</b>
1.1.	Description .....	32
1.2.	Expected use.....	35
1.3.	Therapeutic indications.....	36
1.4.	Contraindications .....	36
1.5.	Precautions for use.....	36
<b>2.</b>	<b>Warnings, symbols.....</b>	<b>36</b>
2.1.	Warnings and safety remarks.....	36
2.2.	Safety linked to the light device .....	37
2.3.	Product labelling.....	38
2.4.	Laser warning sign label .....	39
2.5.	Transport and storage markings .....	39
<b>3.</b>	<b>Technical characteristics .....</b>	<b>39</b>
3.1.	Description .....	39
3.2.	General information .....	40
3.3.	Electrical and optical performances.....	40
3.4.	Physical characteristics of the light source .....	40
3.5.	Ambient conditions .....	41
3.6.	Transport and storage conditions .....	41
3.7.	Inspection of the system .....	42
3.8.	Preparation of the source before treatment .....	42
3.9.	Connection of the optical fibres .....	42
<b>4.</b>	<b>Protocol for use .....</b>	<b>43</b>
4.1.	Preparation of the patient's skin.....	43
4.2.	Installation of the FLUXMEDICARE® system.....	43
4.3.	Putting the intermediate unit in place .....	43
4.4.	Installation of the other accessories FLUXMEDICARE ® .....	45
4.5.	Using the light source .....	46
4.5.1.	Use of the Use of the source .....	47
4.6.	Emergency stop button .....	47
4.6.1.	Display screen.....	47
4.6.2.	Initialisation and choice of language .....	48
4.6.3.	Starting treatment.....	48
4.6.4.	Pausing .....	49
4.6.5.	Stopping at the end of treatment.....	49
4.7.	End of treatment .....	49
<b>5.</b>	<b>Error messages .....</b>	<b>49</b>

<b>6.</b>	<b>Cleaning .....</b>	<b>50</b>
6.1.	General .....	50
6.2.	Source .....	51
6.3.	Cap .....	51
6.4.	The stretch net .....	51
6.5.	Light emitting textiles.....	51
<b>7.</b>	<b>Maintenance .....</b>	<b>52</b>
<b>8.</b>	<b>Recycling .....</b>	<b>52</b>
<b>9.</b>	<b>Electromagnetic compatibility.....</b>	<b>52</b>
<b>10.</b>	<b>Regulation.....</b>	<b>53</b>
<b>11.</b>	<b>Manufacturer details .....</b>	<b>54</b>

# Foreword

---

Thank you for buying the FLUXMEDICARE® system photodynamic therapy system.

Please read this manual carefully before using the FLUXMEDICARE® system and keep it for future reference.  
The content of this manual is based on information updated after compilation.

This manual of FLUXMEDICARE® is available on the product website by following the link:  
<https://www.fluxmedicare.fr/traitemet-pdt/>

## 1. Introduction

---

### 1.1. Description

---

The FLUXMEDICARE® system is a light device intended to be used solely for photodynamic therapy (=PDT) for skin lesions.

The FLUXMEDICARE® system is made up of a red-light generator - called the light source - and a band of luminous textile transmitting the light (light emitting textile).

The source and the light emitting textile are coupled.

The source is connected to the mains electricity supply.

The light emitted has a peak wavelength of 638 nm.

The irradiance received by the patient for the 2.5 hours' treatment is 12 J/cm<sup>2</sup>.

The system is made up of two kits sold separately:

- The **FLUXMEDICARE® No.1 kit** - Ref. LIGHTSOURCE01 - including the light source which is made up of the 3 following elements:
  - A light source
  - A harness (HARN1) for attachment to the torso to hold the intermediate unit. This can be sold separately.
  - A cap (CAP2) which can be sold separately.

Reference	Description	Pictures
LIGHTSOURCE01	1x source 1x harness 1x cap	

The following references are also sold separately:

Reference	Description	Pictures
HARN1	1 x harness	
CAP2	1 x cap	


- The **FLUXMEDICARE® No.2 kit** - Ref. FLUXMEDx - including the light device is made up of the following elements:
  - 12 stretch nets, 15 cm in length (NET57) for single use only,
  - 2 stretch fixing nets specially for the head (NETH) for single use only,
  - 1 batch of 10 pieces of “FLUXMEDICARE® Protect” protective covers (FMP0X) for single use,

- A band of light-emitting textile covered with opaque blue fabric (sheath) and extended by 1 to 3 cables to be connected to the source (TEXTLIGHTx),
- A specific count key to be connected to the source.

Reference	Description	Use	
NET57	6 nets: size No.7 + 6 nets: size No.5	Single use only	
NETH	2 special head stretch nets - size L	Single use only	
<b>TEXTLIGHTX: Where X is a number from 0 to 9 corresponding to the model</b>	1x light emitting textile (the different models are listed on the below table) 1x count key	Re usable	




Reference	Description Light surface	Pictures
TEXTLIGHT1	Bilateral sheath 18 cm x 20 cm and 1x count key	
TEXTLIGHT2	Unilateral sheath 18 cm x 20 cm and 1x count key	
TEXTLIGHT3	Small unilateral sheath 6 cm x 20 cm and 1x count key	
TEXTLIGHT4	Small bilateral sheath 6 cm x 20 cm and 1x count key	

Reference	Description Light surface	Pictures
TEXTLIGHT6	Bilateral sheath for temples 6 cm x 9 cm* and 1x count key	

TEXLIGHT6: schema of light surface

Accessories are available per unit :

Reference	Description	Use	Non contractual pictures
NET57	6 nets: size No.7 + 6 nets: size No.5	Single use only	
NETH	2 special head stretch nets - size L	Single use only	
FMP01	Non woven protective cover for light-emitting textiles TEXTLIGHT1 – Single size	Single use	
FMP02	Non woven protective cover for light-emitting textiles TEXTLIGHT2 – Single size	Single use	
FMP03	Non woven protective cover for light-emitting textiles TEXTLIGHT3 – Single size	Single use	
FMP04	Non woven protective cover for light-emitting textiles TEXTLIGHT4 – Single size	Single use	
FMP06	Non woven protective cover for light-emitting textiles TEXTLIGHT6 – Single size	Single use	

## 1.2. Expected use

The FLUXMEDICARE® system is intended to be used by health professionals practising in doctors' surgeries, hospitals and clinics for PDT (PhotoDynamic Therapy) treatments.

The treatment must be realized only by a medical authorized professional and trained in the PDT.

For the PDT treatment, and before exposure to the light, the area of skin to be treated is covered by a gel or cream containing 5-aminolevulinic acid or an ALA ester or methyl aminolevulinate (not supplied with the FLUXMEDICARE® system) acting as a photosensitiser. These creams are prescription-only.

### **1.3. Therapeutic indications**

---

PDT is used to treat precancerous skin lesions (actinic keratosis) thin and non-pigmented of the face and scalp. (Olsen grade 1 to 2).

The target population is only adults.

### **1.4. Contraindications**

---

- Product contraindicated for pregnant women.
- Product contraindicated for patients with a history of abnormal photosensitivity.
- Patients allergic to the gel or cream containing 5-aminolevulinic acid or an ALA ester or methyl aminolevulinate or the excipients in the cream (not supplied), including ground-nut oil, peanuts, or soya.
- Porphyria
- An immunosuppressed patient for therapeutic reasons, idiopathic or specific disease

### **1.5. Precautions for use**

---

The untreated healthy skin beside the lesion does not need to be protected during exposure.

The FLUXMEDICARE® system is intended to be used in combination with a pharmaceutical cream containing 5-aminolevulinic acid or an ALA ester or methyl aminolevulinate. (not supplied with the FLUXMEDICARE® system)

For contraindications and adverse reactions please refer to the instruction leaflet for these medicines.

The FLUXMEDICARE® system must be used with a transparent occlusive dressing (not supplied with the FLUXMEDICARE® system). This should carry the CE marking.

For contraindications and adverse reactions please refer to the instructions for use of this medical device.

The treatment should only be carried out by an authorised medical professional trained in photodynamic therapy.

## **2. Warnings, symbols**

---

### **2.1. Warnings and safety remarks**

---

This manual contains several warnings describing the banning of or the obligation for certain actions required on the part of the user in order to operate the system in complete safety.

In this manual the warnings are notified by a pictogram  (exclamation mark within a yellow triangle). **Non-observance of these directions for use may lead to serious injury, to the patient or the user!**

In this manual, the information is indicated by a pictogram  (the letter i contained in a blue circle).  
**The following warnings should be read carefully!**

**CAUTION:**

- Do not use the device before reading the user manual carefully and understanding it.
- Do not use electrical power outside the nominal limits specified in the user manual or printed on the back of the device.
- Risk of eye injuries. Check that safety instructions and principles are followed, and there is no possibility of direct light transmission outside the application area.
- Do not look directly at the light emitted.
- Repair and maintenance of the light source is only allowed by MDB TEXINOV personnel or by a distributor authorized by MDB TEXINOV.  
Maintenance of the light source is only allowed by MDB TEXINOV personnel or by MDB TEXINOV approved distributor.
- Prohibition on modification or reparation of the FLUXMEDICARE® system.
- Do not use for indications other than those specified in this document.
- Risk of eye injuries. Do not leave the light fibres illuminated when these are not being applied to the patient.
- Do not connect accessories other than those specified.
- Do not obstruct the air vents or place the appliance on a soft or sloping surface.
- Adhere to a free space or a distance from the wall of at least 30 cm on all sides.
- Do not use an electrical extension lead.
- Only use the count key supplied.
- Do not connect another appliance or count connector (USB type) on the front face of the light source.
- Do not connect this key to a computer or any other electronic appliance.
- Do not place the source on an unstable surface.
- In order to prevent the risk of electric shock, the FLUXMEDICARE® system must be connected to an earthed electric socket.
- The source must not come into contact with a liquid.
- The FLUXMEDICARE® system must not be exposed to high pressures or pointed objects.
- The current socket plug is the disconnecting means and should remain accessible.
- Do not use beyond the duration for use.
- In the treatment room, avoid storing and using flammable materials such as oxygen and some anaesthetics.

## 2.2. Safety linked to the light device

The slightest electronic fault deactivates the light emission.

If the connectors disconnect, the light stops, the source emits a sound signal and the treatment is stopped.

The source complies with the following standards:


- EN NF 60601-1 for electrical safety
- EN NF 60601-1-2 for electromagnetic compatibility
- EN NF 60825 for the safety of laser equipment.

## 2.3. Product labelling

Symbol	Meaning
	Manufacturing date
	European Declaration of Conformity
	Batch number
	Reference
	Medical Device
	Manufacturer
	Type B part applied
	Consult the instructions for use
	Consult the precautions for use
	Laser class 1
	Temperature limit
	Compliant with directives relating to Waste Electrical and Electronic Equipment (WEEE) and to the limitation of certain hazardous substances in electrical and electronic equipment (RoHS). This symbol affixed to the product or its packaging indicates that the product cannot be discarded in ordinary household waste.
	Expiry date
	ON button (turning on)
	OFF button (turning off)
	Single use: do not re-use the net

## 2.4. Laser warning sign label

The following label is located on the rear part of the laser (original black on yellow). It shows the laser class.

Class1 laser emission

## 2.5. Transport and storage markings

Symbol	Meaning
	The arrows point upwards, store and transport the parcel vertically
	Acceptable temperature limits (-25°C to 50°C)
	Fragile, avoid knocking and jolting
	Store in a dry place
	Acceptable range of ambient pressure
	Acceptable relative humidity range

## 3. Technical characteristics

### 3.1. Description

The power emitted does not need to be calibrated by the user because it is controlled by the device itself.

### 3.2. General information

---

Parameters	Unit	Values	Comments
<b>Wavelength</b>	nm	638	-
<b>Irradiance</b>	mW/cm <sup>2</sup>	1,3	On the light bands: > 1 mW/cm <sup>2</sup>
<b>Treatment duration</b>	min	150	± 1 continuously

Description	Classification
<b>Laser type</b>	Diode laser
<b>Laser class</b>	1
<b>Protection against electric shocks</b>	Class I
<b>Protection against harmful penetration of water or particles</b>	IP X0
<b>Part applied: band of luminous textile</b>	Type B
<b>Cooling system</b>	PELTIER


### 3.3. Electrical and optical performances

---

Parameter	Range	Unit
<b>Input voltage</b>	220-240	VAC
<b>Input frequency</b>	50-60	Hz
<b>Primary power</b>	250	W
<b>Input current</b>	3	A
<b>Output voltage</b>	12	V
<b>Output current</b>	17.5	A

### 3.4. Physical characteristics of the light source

---

### 3.5. Ambient conditions

The system is designed to be used in a standard hospital or clinic environment, with the following ambient conditions:

Parameter	Unit	Min	Max	Comment
<b>Temperature for use</b>	°C	+15	+30	
<b>Humidity for use</b>	% RH	10	85	Environment for use without condensation required
<b>Atmospheric pressure</b>	hPa	700	1060	

### 3.6. Transport and storage conditions

Parameter	Unit	Min	Max
<b>Transport and storage temperature</b>	°C	-25	+50
<b>Transport and storage humidity</b>	% RH	10	85
<b>Atmospheric pressure</b>	hPa	700	1060

After transport or storage, the user must ensure the device is installed in its normal conditions of use (see §3.5. Ambient conditions).

The system must not be stored near a heat source or in direct sunlight.

The light emitting bands must be stored flat; the optical fibres must not be folded or crushed.

Avoid exposure to very dusty environments. Avoid storage and use of flammable materials such as oxygen and some anaesthetics in the treatment room.

The intermediate unit can be maintained on the main unit thanks hook-and-loop fastening strips.

Caps are placed on the connector ends of the intermediate unit.

### **3.7. Inspection of the system**

**i** After unpacking the FLUXMEDICARE® device, inspect it and check it has not suffered any damage during transport.

Before use, make sure the device has not been damaged, that the various parts are intact, e.g.

- The various parts are not broken
- The seams on the cap are in good condition
- The optical fibres are not damaged

If there is a damage, contact MDB Texinov before use.

### **3.8. Preparation of the source before treatment**

#### **Environment**

**i** The equipment must be used on the premises of a health professional. Do not use in conditions where the ambient temperature exceeds 30°C or where there is condensation.

#### **Installation**

The FLUXMEDICARE® source from MDB Texinov must be installed on a stable flat surface. Particular attention must be paid to checking that the ventilation outlets are not obstructed and that sufficient space between the front and rear panels and the adjacent surface is maintained (at least 10 cm), in order to prevent any obstruction to the air flows that could cause internal overheating of the system.

**i** These caps must be removed before the connection.

### **3.9. Connection of the optical fibres**

The connection between the MDB Texinov source and the optical fibres must be made before starting the system.

All connectors must be connected to the intermediate unit.

## 4. Protocol for use

---



### 4.1. Preparation of the patient's skin

Before applying the photosensitising cream, the surface of the lesions should be prepared to remove the scales and scabs and to roughen the surface of the lesions.


Using a spatula apply a layer of photosensitising cream (about 1 mm thick) on a radius of 5 to 10 mm of healthy skin around the lesions. Cover the area to be treated with a transparent occlusive dressing.

### 4.2. Installation of the FLUXMEDICARE® system

1. Place the light emitting textile (sheath) on the patient's head. The light emitting textile is put in place with the light emitting face above the lesions to be treated and on the transparent dressing.
2. Position the stretch net on the navy-blue side of the sheath so that the light emitting bands are fitted on the areas to be treated and to prevent the light from escaping.
3. Put the cap on the net. The cap can be put on in any direction. The system's various layers must remain positioned as they are for a period of 30 minutes before the start of the light emission. After 30 minutes the adjustment of the cap is checked before starting the light emission.

### 4.3. Putting the intermediate unit in place

### **Putting the harness and the intermediate unit in place**

1. Place the light source on a stable surface, e.g. a table beside the person to be treated.

2. Put the chest harness in place.

The attaching harness is completely adjustable. The adjustment is made using the 5 adjustment buckles.

- Open the attaching harness using the clip by pressing it.

- Put the harness on like a waistcoat.

- 

- Close the harness again with the clip.

- Tighten or loosen the buckles on the chest and waist to adjust it.

- Unscrew the locking screw.

- Insert the intermediate unit.

- Retighten the locking screw.
- Check the intermediate unit is correctly attached.

3. Connect the three connectors to the intermediate unit. These should be above the harness.

Push the metal end (connector) of the light emitting textile into the intermediate unit.

4. Connect the FLUXMEDICARE® count key supplied with the light emitting textile.

This must remain connected throughout the treatment.

A connection is provided on the face of the light source (COUNT KEY) .


#### 4.4. Installation of the other accessories FLUXMEDICARE ®

1. Place the light emitting textile on the head or the face of the patient.  
The transparent wound dressing is positioned on the to the lesions to be treated.  
the light emitting bands are put with the light emitting face above the lesions to be treated and on the transparent wound dressing.
2. Position the stretch net on the navy-blue side of the sheath so that the light emitting bands are fitted on the areas to be treated and to prevent the light from escaping.
-

Reference	Photography	Reference	Photography
<b>TEXLIGHT2 or TEXLIGHT3 Unilateral sheath</b>			
<b>TEXLIGHT6 bilateral sheath for temples</b>			

#### 4.5. Using the light source

Diagram 2: Front face of the light source

The front face of the light source (see Diagram 1) is made up of:

- A display screen (user interface)
- A START/PAUSE button
- A STOP button
- A connection for the count key (COUNT KEY)
- A connection for the reinforced cable which connects the source to the intermediate unit
- A yellow emergency stop button, to be used to stop the source in an emergency (EMERGENCY STOP)
- An ON (I)/OFF (O) start button
- A green indicator light when the source is on (ON/OFF)
- A red indicator light which comes on if the source is overheating, or if an emergency stop is triggered (EMERGENCY STOP/OVERHEATING)
- A blue indicator light which comes on when the source is powered (POWER).

#### 4.5.1. Use of the Use of the source

1. Connect the light source to the mains electricity.  
A blue indicator light comes on when the appliance is powered (ON/OFF light).
  
2. Turn on the source.  
Press the rocker button ON (I)/OFF (O).  
The display screen comes on.

#### 4.6. Emergency stop button

The emergency stop switch immediately stops the emission of the light beam. It is positioned on the front face of the light source by the EMERGENCY STOP label (yellow/red button).

##### Activate the emergency stop

- Press the EMERGENCY STOP button (yellow).
  - The emission of the laser light beam is stopped immediately.
  - The display screen goes off.
  - The appliance is not then turned off from the mains, the POWER light remains blue.
  - The "EMERGENCY STOP/OVERHEATING" indicator light (red) comes on.
  - The "ON/OFF" light (green) remains lit.

##### Disarming the emergency stop

- Turn the emergency stop switch to the right until the stopping device snaps in place and the red switch is triggered forwards.
  - It is still possible to use the appliance.
  - The "EMERGENCY STOP/OVERHEATING" indicator light goes off.
  - The display screen comes on.


#### 4.6.1. Display screen

The display screen is on the front panel. The verification of parameters and the starting of treatment take place through this interface via the buttons located on the front of the appliance.

The screen displays all the important parameters related to the treatment and the device functions.

The light emission may be paused and restarted via the button (START/PAUSE); the status is shown on the user interface.


#### 4.6.2. Initialisation and choice of language


- During the display of the initialisation screen, the system carries out a self-diagnosis and configures the source parameters. The software version is specified on the screen when the source starts up.

During this display it is possible to choose the software language (German, English, French):


- Press the START/PAUSE and STOP buttons simultaneously
- Scroll down the languages by pressing the START/PAUSE and STOP buttons simultaneously
- Confirm the language chosen by a long press (>5 s) on the START/PAUSE button
- Cancel the language chosen by a long press (<5 s) on the STOP button
- Once the language is chosen it will remain in the memory. The previous step is no longer necessary.


#### 4.6.3. Starting treatment

- After initialisation, the warm-up of the lasers is carried out.

- The connected lasers are shown ON. Check that all the selected lasers are properly connected.

- Press the START/PAUSE button to continue.

- The screen displays the number of treatments remaining.

- Press the START/PAUSE button to continue.

- Press the START/PAUSE button to start the treatment.

- The treatment starts.

The remaining time is shown on the screen (in h: min: s).

#### 4.6.4. Pausing

---

- When the lights emit the screen displays the remaining treatment time.




- The treatment can be stopped for a short while by pressing the START/PAUSE button. A beep sounds when pressed.
- After pressing the START/PAUSE button again, the treatment re-starts. This pause is limited to 15 min.
- If this time is exceeded the treatment stops.

#### 4.6.5. Stopping at the end of treatment

---


- When the lights emit the screen displays the remaining treatment time.

- The treatment can be stopped by pressing the STOP button. A beep sounds when pressed. The light emission stops.
- Press the STOP button a second time to stop the treatment. A sound signal is heard until the source is turned off by the user with the ON (I)/OFF (O) button.
- Then the screen indicates the end of treatment and then that the treatment has been cancelled.

#### 4.7. End of treatment

---

- The screen shows the end of treatment.

- At the end of treatment, the light emission stops. A sound signal is heard until the source is turned over by the user.
- Use the ON (I)/OFF (O) button to turn the source off.
- Remove the cap, the stretch net, the light emitting textile and then the transparent dressing. The transparent dressing and the stretch net must be discarded in the appropriate bin.
- Clean the area of the skin with saline solution.

### 5. Error messages

---

NB: X is a number from 1 to 3 corresponding to the 3 lasers of the light source

Error	Description	Action required
1.	Temperature error on X laser Previous treatment error	It appears on starting when an error occurred during the previous treatment. Delete the message by pressing STOP; keep it by pressing START/PAUSE. An operating diagnostic

Error	Description	Action required
	START: Memorise STOP: Delete	will then be performed at the next start-up. If the problem persists contact Technical Support.
2.	Irradiance error on laser X Previous treatment error START: Memorise STOP: Delete	It appears on starting when an error occurred during the previous treatment. Delete the message by pressing STOP; keep it by pressing START/PAUSE. No diagnostic performed. If the problem persists contact Technical Support.
4.	Light emitting textile no longer usable - it must be replaced.	Contact MDB Texinov to obtain new light emitting textiles.
5.	ERROR: Please choose at least 1 laser START: <—	Connect at least one connector of the light emitting textile to the source and press the START/PAUSE button.
6.	! LASER X ERROR! Connect laser and press START.	Connect the connector of the light emitting textile (1, 2 or 3 of the intermediate unit) to the source and press the START/PAUSE button. While the laser is not reconnected a sound, signal is heard.
7.	! LASER ERROR! Connect lasers and press START.	Connect the connector of the light emitting textile to the intermediate unit and press the START/PAUSE button. While the laser is not reconnected a sound, signal is heard.
8.	/!\ ERROR /!\ Thermal anomaly detected on LASER X	Stop the source and wait until it cools. Check that the air openings are not obstructed and that the ambient temperature does not exceed the nominal values. If the problem persists contact Technical Support. When starting next the appliance will display an error. A diagnostic will be performed.
9.	/!\ ERROR /!\ Irradiance anomaly detected on LASER X	Stop the source and wait until it cools. If the problem persists contact Technical Support. When starting next the appliance will display an error A diagnostic will be performed.
10.	! CONNECTION ERROR! Connect the memory And Press START	The count key is not connected. Connect the count key and press START/PAUSE.
11.	! MEM ERROR! Serial No. not valid Connect and press START	The count key is not the one that was connected at the start of treatment. Reconnect the count key used initially and press START/PAUSE.

## 6. Cleaning

---

### 6.1. General

---

- CAUTION:**
- Any cleaning must be carried out when the source is not connected to the mains electricity.

- The light emitting textile is disconnected from the source.
- NEVER IMMERSE a photodynamic therapy component in a liquid.
- The device must not be sterilised by autoclave or gas.
- Do not bleach.
- Do not dry in a machine.
- Do not iron.
- Do not dry clean.
- Wash only with soapy water.

## 6.2. Source

---

- After turning off the source wait 5 minutes before cleaning.
- Disconnect the power supply lead for the source before cleaning (risk of electric shock).
- Wet a soft, clean cloth with water. Clean the outside of the source and the power supply lead with the cloth. Make sure that no liquid enters the source.
- Dampen a second soft, clean cloth with water. Wipe the source.
- Before use make sure that the source is completely dry. Check that the electric lead is not damaged.

## 6.3. Cap

---

The cap can be washed in the machine on the delicate programme at 40°C without fabric softener.

## 6.4. The stretch net

---

After use discard the single use net.

## 6.5. Light emitting textiles


---

In case of dirty light strips please follow the instructions below:

- Clean with soapy water
- let the material air dry. DO NOT DRY WITH ARTIFICIAL HEAT.
- Wipe the light strips with a dry cloth.

Before each use, check the condition of the light strips. Particular attention should be paid to:

- The state of cleanliness: If, despite the cleaning described above, the light strips are stained or dirty, they must be discarded.
- The condition of the fibres on the luminous part. Light strips must be discarded if they have more than 10 light points (see photo below)

## 7. Maintenance

---

The source must be checked after each use for obvious visible damage, damaged cables, cracks, worn label.

The source is guaranteed for one year against material or manufacturing defects, from the delivery date. During this period, the defective equipment will be replaced free of charge, unless the fault is due to abuse or incorrect use.

The source must be checked annually by MDB Texinov® personnel.

This control can also be carried out by a distributor approved by MDB Texinov®.

The service life of the light source is 5000 hours.

MDB Texinov will only take responsibility for the safety, reliability and perfect operation of the appliance if:

- All the adjustments as well as modifications or repairs have been carried out by MDB Texinov.
- The appliance is used in compliance with the instructions for use supplied.

### **Important indications**

Only MDB Texinov technicians are authorised to repair FLUXMEDICARE® products.

The units must not be opened.

Any serious incident in connection with the device should be notified to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established.

## 8. Recycling

---

Electrical appliances must not be discarded with household waste.

In France, it must be taken to an appropriate collection centre for recycling of electrical and electronic equipment.

To ensure correct treatment of product waste, refer to the usual recycling network, in compliance with local requirements for electrical equipment.

## 9. Electromagnetic compatibility

---

- The use of this appliance beside other appliances or stacked with them should be avoided. This could cause a malfunction. If this use is necessary this appliance, and the others, should be monitored to ensure they are functioning normally.
- The use of accessories, transducers and cables other than those specified or supplied, may cause an increase in electromagnetic emissions or a decrease in the immunity of this machine and cause inappropriate operation.

- RF portable communications devices (including peripherals such as antenna cables and external antennae) should not be used closer than 30 cm (12 inches) from any part of the APPLIANCE, including cables specified by the manufacturer. Otherwise the performance of these devices could be affected.

Electromagnetic compatibility (EMC) and electrostatic discharge (ESD) tests have been performed in accordance with the following standards:

Phenomenon	Standard base used	Result
<b>Radiated emission</b>	CISPR 11+A1: 2010	Meets the requirements of Group 1 Class B
<b>Emission by conduction</b>	CISPR 11 Group 1 Class B	Meets the requirements of Group 1 Class B
<b>Harmonic current emission 50 Hz</b>	EN/IEC 61000-3-2+A1+A2: 2009	Meets the requirements of Class A
<b>Voltage fluctuation</b>	EN/IEC 61000-3-3: 2013	Meets the requirements
<b>Electrostatic discharges</b>	IEC 61000-4-2: 2008	Meets the requirements of Criterion A
<b>Radiated electrostatic fields</b>	EN 61000-4-2: 2008	Meets the requirements of Criterion A
<b>Fast transient electrical disturbance/burst</b>	EN/IEC 61000-4-4-3+ A1+A2: 2010	Meets the requirements of Criterion A
<b>Shock wave</b>	IEC 61000-4-5: 2005	Meets the requirements of Criterion A
<b>RF conducted disturbances</b>	EN 61000-4-6: 2013	Meets the requirements of Criterion A
<b>Magnetic field 50 Hz</b>	IEC 61000-4-8: 2009	Meets the requirements of Criterion A
<b>Supply voltage interruption and dip</b>	EN/IEC 61000-4-11 :2004	Meets the requirements of Criterion A

## 10. Regulation

The laser source:

- is a Class IIa appliance compliant with the regulations on medical products
- is a laser Class 2M laser device (in accordance with the standard, NF EN 60825-1)
- has CE marking in accordance with Directive 93/42/EEC.

The FLUXMEDICARE® V2 device is a CE marked Class IIa medical device in accordance with Directive 93/42/EEC obtained in 2018 by the BSI organisation (CE2797).

The accessories are CE marked Class I medical devices in accordance with Regulation (EU) 2017/745.

The following standards are applicable.

Reference	Title
<b>NF EN ISO 13485 version: September 2016</b>	Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes

Reference	Title
<b>Directive 93/42/EEC</b>	Directive 93/42/EEC from the Council meeting of 14 June 1993 relating to medical devices, amended by Directive 2007/47/EC
<b>Regulation (EU) 2017/745</b>	New European regulation of 05 April 2017
<b>NF EN 60825-1 version: October 2014</b>	Safety of Laser Products - Part 1: Equipment classification, requirements and user's guide
<b>NF EN ISO 15223-1 version: January 2017</b>	Symbols to be used with medical device labels, labelling and information to be supplied
<b>NF EN 60601-1/A1 version: March 2014</b>	Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance
<b>CEI 60601-1-2:2014 version: February 2014</b>	Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
<b>NF EN 60601-1-6 version: July 2010</b>	Medical Electrical Equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
<b>NF EN 62304 version: October 2006</b>	Medical device software - Software life-cycle processes
<b>NF EN 62304/A1 version: April 2018</b>	Medical device software - Software life-cycle processes Amendment A1
<b>NF EN 62366-1 version: December 2015</b>	Medical devices - Part 1: Application of usability engineering to medical devices
<b>NF EN ISO 14971 version: December 2019</b>	Medical devices - Application of risk management to medical devices

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